

## Human anti- Staph Enterotoxin B mAb (628)

**Catalog Number:** S200-008

**Lot:** 2204006

**Concentration:** 100 µg is supplied in PBS at 2.6 mg/mL.

**Description:** *Staphylococcus aureus* (*S. aureus*) is a Gram-positive bacterium that colonizes humans and domestic animals. *S. aureus* produces over 20 different enterotoxins, of which Enterotoxin A and B (SEB and SEA) are the most well characterized and also regarded as superantigens due to their ability to bind to MHC class II alpha and beta chains. SEA is the most common toxin associated with staph food poisoning, and SEB has been studied as its potential as an inhaled bioweapon. Eliteimmune's SEB targeting mAb 628 can be used in toxin neutralization assays (TNA) and is also a useful research tool for ELISAs.

**Product Type:** Monoclonal Antibody IgG1 subtype

**Strain:** *Staphylococcus aureus*

**Synonyms:** Enterotoxin B, SEB, Staphylococcal enterotoxin B, superantigen B, superantigen

**Expression Host:** CHO

**Tag:** None

**Purity:** Antibody is purified using immobilized protein G.

**Applications:** ELISA and toxin neutralization assays (TNA). Specific conditions should be optimized by end user.

**Formulation:** PBS

**Storage:** 4°C for immediate use, -80°C for long-term storage. Avoid repeated freeze thaw cycles.

**Shipped:** Dry Ice

**Intended Use:** For research use only. Not intended for human, therapeutic, or diagnostic use. Buyer may not modify, sell, or transfer product for commercial use without written permission from Eliteimmune.

### Recommended Dilutions:

**ELISA:** Assay-dependent dilution. The antibody can be used for capture of SEB in an indirect ELISA assay.

**WB:** Not recommended for Western Blot because the antibody recognizes conformational epitope.

**Neutralization:** Assay-dependent dilution

**Cross Reactivity:** This antibody does not appear to cross react with *S. aureus* antigens Hla, LukF, LukS, or SEA at the tested concentrations.

### ELISA Data:

ng/mL	OD 650 nm				
	Hla	LukF	LukS	SEA	SEB
100.000	0.096	0.084	0.094		3.703
10.000	0.092	0.080	0.076	0.150	3.404
1.000	0.086	0.083	0.088	0.083	3.456
0.100	0.090	0.074	0.081	0.064	3.690
0.010	0.083	0.079	0.079	0.072	3.402
0.001	0.090	0.083	0.078	0.093	1.674
0.000	0.088	0.071	0.119	0.107	0.405
0.000	0.092	0.074	0.081	0.100	0.138

*S. aureus* antigens were diluted to 2 µg/mL in PBS for plate coating. Anti-SEB antibody 628 was serially diluted ten-fold from 100 µg/mL and incubated on the coated plates. Washed plates were detected with anti-human IgG-HRP conjugate and TMB substrate. OD650 is reported above.